



## INSTRUCTIONS FOR TESTING, CONSULTATIONS, AND SUBMITTING DIAGNOSTIC SAMPLES

√ Test	Test	Instructions	Cost US\$
	THYROID ANTIBODY PROFILE (D8T) Test includes; T3, T4, free T3, free T4, T3 autoantibodies and T4 autoantibodies.	If the patient is already on thyroid medication the blood should be drawn 4-6 hours after giving thyroid medication. Collect 3 to 5 mL of blood. Add to redtop tube or serum separator tube (SST). Allow to clot and centrifuge to separate serum. <b>DO NOT SEND UNSEPARATED BLOOD.</b> Transfer serum to empty redtop or plastic tube. 1.0 mL serum is <b>REQUIRED</b> for the Thyroid Antibody Profile. Test includes T3, T4, free T3, free T4, T3 autoantibodies and T4 autoantibodies.	\$65.00
	THYROID ANTIBODY PROFILE + THYROGLOBULIN AUTOANTIBODY (D8T + TgAA)	Same as above but send extra 0.5 mL separated serum	\$85.00
	Thyroid Ab Profile, TgAA and OFA REGISTRY PANEL	See above for guidelines except – 2.0 mL serum is <b>REQUIRED</b> to complete all testing. <b>For OFA, include completed &amp; signed OFA form, and submit separate check for \$15.00 made out to OFA.</b>	\$ 85.00
	PROFILE 7200 includes CBC, Differential & Chemistries & D8T	Send 1 mL whole blood in lavender-top tube (EDTA) <b>AND</b> 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack when shipping, but do not freeze sample. Without the CBC, deduct \$6.	\$87.00
	TgAA (thyroglobulin autoantibody), Add -On to D8T	<u>Note:</u> Good for genetic screening. Send extra 0.5 mL separated serum. Allow <b>3 - 5 days testing time.</b>	\$27.00
	T4 and freeT4 <u>only</u> (44)	Send 0.5 mL separated serum	\$51.00
	CBC, Differential & Chemistries (D2)	Send 1 mL whole blood in lavender-top tube (EDTA) <b>AND</b> 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack for shipping, but do not freeze. Minus CBC, deduct \$6.	\$52.00
	DISTEMPER & PARVO TITERS	Send 1.0 mL separated serum.	\$35.00
	RABIES VACCINE TITER	Send 1.5 mL serum in red top or plastic tube. <u>Note:</u> Allow <b>14 - 21 days testing time.</b>	\$75.00
	VON WILLEBRAND TEST wWD	Take 1 mL separated plasma from a citrate tube (BTT, light blue) and transfer plasma to empty redtop or plastic tube. Include ice pack when shipping, but do not freeze sample. Clearly label the tube as "PLASMA". <u>Note:</u> Allow <b>4-7 days testing time.</b>	\$55.00
	PHENOBARBITAL PROFILE	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. <u>Please use red top tube, NOT an SST tube.</u>	\$45.00
	BROMIDE PROFILE	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. <u>Please use red top tube, NOT an SST tube.</u>	\$53.00
	<b>SERUM BILE ACIDS PRE &amp; POST</b>	1st tube is fasting, 2 nd tube is 2-4 hours post-fatty meal. 1.0 mL separated serum in <b>each</b> tube. Label tubes clearly 1 st & 2 nd.	\$44.00
	<b>URINE BILE ACID</b>	Freshly collected urine, 3 mL.	\$43.00
	PANLEUKOPENIA VACCINE TITER – For <b>Felines</b>	Send 1.0 mL separated serum.	\$35.00
	CONSULTATION REQUEST only	<b>(Attach all lab and/or medical information)</b>	\$ 25.00

- **LABEL TUBES CLEARLY WITH PETS NAME AND OWNERS LAST NAME**
- Ship samples in a **sturdy cardboard box** with plenty of padding.
- **Do not** freeze samples. Blue & lavender top tubes require an ice pack with an overnight courier.
- Send by US Priority Mail, US Express Mail, Fed-Ex, UPS or other courier.
- Samples received in our office Monday – Saturday.
- **HEMOPET** will process and submit samples to Antech Diagnostics for analysis. **Dr. Dodds** will follow with interpretations and recommendations.
- Send check, money order or credit card info along with samples to:

**HEMOPET/HEMOLIFE**  
11330 Markon Drive  
Garden Grove CA 92841

Any questions, please call our office – 714-891-2022 Monday to Friday, 8:30 AM – 5:30 PM.

Download additional forms from <http://www.hemopet.org/services.html>